

PT. KRISAN KALINDO UTAMA

Distributor Alat Kesehatan & Laboratorium Terpercaya

Email: info@krisankalindoutama.com | Telp: +62-21-826-58-722

BIOBASE ELECTROLYTE ANALYZER BKE-B



SKU : BE-02

Kategori : Electrolyte Analyzer

DESKRIPSI PRODUK

The product is used to analyze the concentration of electrolyte ions potassium (K⁺), sodium (Na⁺), chloride (Cl⁻) in Human serum, plasma, and whole blood samples. The structure of the product is composed of electrode module, testing module, pipeline module, circuit module and data output module. BKE-B has a sample tray.



Scan QR Code atau kunjungi link berikut untuk melihat detail di website:

<https://mail.krisankalindoutama.com/products/biobase-electrolyte-analyzer-bke-b>