



## H1650K Benchtop High Speed Centrifuge

H1650K Benchtop High Speed Centrifuge provides multiple protection against over-speed, door lid and misoperation, etc. 7 types of rotors are available for users to choose from, with powerful functions and safe operation, which can meet various centrifugal experiments.



### Product advantages








- ◆ 3.5 inch HD touch screen, simple interface, setting parameters and running parameters can be displayed simultaneously.
- ◆ Directly input the setting value for fast and convenient operation without adjust the parameters back and forth for many times. (compared to traditional keypad and knob control products).
- ◆ Multiple user-defined programs can be saved.
- ◆ Microprocessor control, DC brushless motor for stable and low noise operation.
- ◆ Parameters can be set freely during a run.
- ◆ door lid interlock to improve the safety performance.
- ◆ Compact with great performance.
- ◆ 10 Acceleration and 11 Deceleration rates, to achieve best centrifugation result.



### Technical parameter

Order NO	0203191000
Max Speed	16500r/min
Max RCF	18360xg
Max Capacity	12×5ml
Speed Accuracy	±30 r/min
Time Setting Range	1s ~ 99min 59sec
Noise	≤ 65dB (A)
Power Supply	AC 220±22V 50Hz
Power Consumption	250W
Dimension (WxDxH)	320×370×235(mm)
Packing Dimension (WxDxH)	530×410×290(mm)
Net Weight	11kg

 Rotor parameters

<p>No:1</p>  <p>Rotor Type: Angle rotor Max Capacity: 12×1.5ml Max Speed: 16500r/min Max RCF: 18360xg Order NO: 0301809100 Adapter: 0.5ml, 0.2ml</p>	<p>No:2</p>  <p>Rotor Type: Angle rotor Max Speed: 16500r/min Max RCF: 18360xg Max Capacity: 18×0.5ml Order NO: 0301810500</p>
<p>No:3</p>  <p>Rotor Type: Angle rotor Max Speed: 12000r/min Max RCF: 10142xg Max Capacity: 12×5ml Order NO: 0301809200 Adapter: 0.5ml, 1.5ml</p>	<p>No:4</p>  <p>Rotor Type: Angle rotor Max Speed: 13000r/min Max RCF: 16170xg Max Capacity: 24×1.5ml Order NO: 0301809300 Adapter: 0.5ml, 0.2ml</p>
<p>No:5</p>  <p>Rotor Type: Angle rotor Max Speed: 15000r/min Max RCF: 13800xg Max Capacity: 3×8×0.2mlPCR tube Order NO: 0301809400</p>	<p>No:6</p>  <p>Rotor Type: Angle rotor Max Speed: 12000r/min Max RCF: 14800xg Max Capacity: 12capillaries Order NO: 0301856400</p>
<p>No:7</p>  <p>Rotor Type: Angle rotor Max Speed: 12000r/min Max RCF: 14800xg Max Capacity: 24capillaries Order NO: 0301856300</p>	<p>/</p>



**PT. KRISANKALINDO UTAMA**  
 Jl. Baru Underpass No.8. 001/001, Duren Jaya, Bekasi Timur

Telp : 021-82658722  
 WA Admin : 0813 1986 2036  
 E-mail : [krisankalindo@gmail.com](mailto:krisankalindo@gmail.com)