

L500 Tabletop Low Speed Centrifuge

L500 Tabletop Low Speed Centrifuge is newly designed for the clinical research based on our years of experience. Stainless steel rotors with variety of capacity for convenient operation.



Features

- ◆ Microprocessor control and DC brushless motor.
- ◆ Touch panel and LCD display.
- ◆ More larger space chamber compared with its old version;
- ◆ RCF value can be calculated automatically.
- ◆ Special damping structure to reduce vibration.
- ◆ Electrical door interlock to control the door operation for safety guarantee; Door cannot be opened during working procedure.
- ◆ Variety of stainless steel brackets for diversified choice.

Technical

Model	L500	Time Setting Range	1min to 99min
Order No.	209000	Noise	≤65dB(A)
Max. Speed	5000rpm	Power Supply	AC220V±22V 50/60Hz 5A
Max. RCF	4030xg	Total Power	250W
Max. Capacity	16x15ml	Dimension (W x D x H)	380*500*300mm
Speed Accuracy	±30rpm	Package Dimension(W x D x H)	460*660*350mm
Net Weight	19kg		

L500 Rotor Options

 <p>No1</p> <p>Order No.: 808350 Max. Speed: 5000rpm Max. RCF: 4030xg</p> <p>Rotor Capacity: 24x10ml Adaptor: 808308 1.5ml</p>	 <p>No2</p> <p>Order No.: 855201 Max. Speed: 4200rpm Max. RCF: 2760xg</p>	 <p>Name: Bracket Order No.: 855210 Capacity: 4x50ml Adaptor: 854708 15ml falcon</p>	 <p>Name: Bracket Order No.: 855220 Capacity: 4x4x15ml Adaptor: 854806(1.5ml) 854808(5ml)</p>	 <p>Name: Bracket Order No.: 855230 Capacity: 4x6x5/10ml Adaptor: 854807(1.5ml) 854813(0.5ml)</p>
--	---	---	--	--



PT. KRISANKALINDO UTAMA
Jl. Baru Underpass No.8. 001/001, Duren Jaya, Bekasi Timur

Telp : 021-82658722
WA Admin : 0813 1986 2036
E-mail : krisankalindo@gmail.com