

PT. KRISAN KALINDO UTAMA

Distributor Alat Kesehatan & Laboratorium Terpercaya

Email: info@krisankalindoutama.com | Telp: +62-21-826-58-722

FLUORECARE D-DIMER



SKU : MF-34

Kategori : Immunology Analyzer

DESKRIPSI PRODUK

The fluorecare® D-dimer assay is a fluorescence immunochromatographic test that quantitatively measures D-dimer levels in human plasma. The fluorescence signal is proportional to the D-dimer concentration and is automatically analyzed by the instrument, with results reported in mg/L.

Presentation-style version:

Fluorecare® D-dimer is a rapid fluorescence immunoassay that quantitatively detects D-dimer in human plasma. The analyzer measures fluorescence intensity and converts it into a D-dimer concentration (mg/L).



Scan QR Code atau kunjungi link berikut untuk melihat detail di website:

<https://mail.krisankalindoutama.com/products/fluorecare-d-dimer>