

# PT. KRISAN KALINDO UTAMA

Distributor Alat Kesehatan & Laboratorium Terpercaya

Email: info@krisankalindoutama.com | Telp: +62-21-826-58-722

## SARTORIUS BSA 2245 ANALYTICAL BALANCE



**SKU : TS-01**

**Kategori : Uncategorized**

### DESKRIPSI PRODUK

The Sartorius BSA 2245 delivers high precision with a maximum capacity of 220 grams and a resolution of 0.1 mg (0.0001 g). This balance is ideal for chemical, pharmaceutical, and research laboratories that require extreme accuracy in weighing.

The automatic internal calibration feature (isoCAL) ensures consistent results without the need for manual calibration. Its compact design, featuring an easy-to-read LCD screen, makes it user-friendly for daily operations.

Key Features:

- Capacity & Resolution: 220 g capacity, 0.1 mg resolution
- Calibration: Automatic internal calibration (isoCAL)

- Design: Ergonomic and space-saving design
  - Display: High-contrast LCD screen
  - Connectivity: RS-232 port for data connection
- 



Scan QR Code atau kunjungi link berikut untuk melihat detail di website:

<https://mail.krisankalindoutama.com/products/sartorius-bsa-224s-analytical-balance>